

Standard Primers...

We have a stock of standard primers according to the below standard primer list.
The standard primer list contains the most popular standard primers.

- Deprotected & purified (Reverse Phase Column)
- Recommended for methods like PCR, sequencing and hybridisation
- Ready-to-use for Sanger Sequencing
- Standard delivered as 50µM aqueous solution
- 5 nmol of highly pure primer in 100 µl H₂O
- Certificate of Analysis (CoA) follows
- 100% Quality Controlled (QC) by MassSpec
- To be kept frozen at -18 C
- Avoid multiple freeze-thaw cycles

| Product No. | Name | Sequence (5'-3') | Price (DKK) |
|--------------------|------------------|---------------------------|--------------------|
| TAGC1001 | 16S-1392r | ACGGGCGGTGTGTGTRC | 50 |
| TAGC1002 | 16S-27f | AGAGTTTGATCMTGGCTCAG | 50 |
| TAGC1003 | 3'-Gal 4 | GAAGTTGCGGGGTTTTTC | 50 |
| TAGC1004 | 35S-A | AAGGGTCTTGCGAAGGATAG | 50 |
| TAGC1005 | 35S-B | AGTGGAAAAGGAAGGTGGCT | 50 |
| TAGC1006 | BAC1 | GCTCGTATGTTGTGTGGAATTG | 50 |
| TAGC1007 | BAC2 | CCCAGTCACGACGTTGTAAAAC | 50 |
| TAGC1008 | BGH-Reverse | TAGAAGGCACAGTCGAGG | 50 |
| TAGC1009 | CMV-F | CGCAAATGGGCGGTAGGCGTG | 50 |
| TAGC1010 | EBV-RP | GTGGTTTGTCCAAACTCATC | 50 |
| TAGC1011 | EGF1 | CCAAGCTATTTAGGTGACAC | 50 |
| TAGC1012 | EGF2 | GGGGATGTGCTGCAAGG | 50 |
| TAGC1013 | GFPN1-5 | GAGCTCAAGCTTCGAATTC | 50 |
| TAGC1014 | GLprimer1 | TGTATCTTATGGTACTGTAAGT | 50 |
| TAGC1015 | KAN2-FP | ACCTACAACAAAGCTCTCATCAACC | 50 |
| TAGC1016 | KAN2-RP | GCAATGTAACATCAGAGATTTTGAG | 50 |
| TAGC1017 | KS | TCGAGGTCGACGGTATC | 50 |
| TAGC1018 | LPP1 | ATACGACTCACTATAGGGAGAC | 50 |
| TAGC1019 | LPP2 | ACAAATTGGACTAATCGATGGC | 50 |
| TAGC1020 | M13-FP | TGTAAAACGACGGCCAGT | 50 |
| TAGC1021 | M13-RP | CAGGAAACAGCTATGACC | 50 |
| TAGC1022 | pACT2-FP | GATGATGAAGATACCCAC | 50 |
| TAGC1023 | pACT2-RP | CAGTTGAAGTGAAGTTGC | 50 |
| TAGC1024 | pBacPAC-RP | GTCTGTAAATCAACAACGC | 50 |
| TAGC1025 | pBAD-FP | ATGCCATAGCATTATCC | 50 |
| TAGC1026 | pBakPAC-FP | TAAAATGATAACCATCTCGC | 50 |
| TAGC1027 | pBR1 | CGAAAAGTGCCACCTGAC | 50 |
| TAGC1028 | pBR2 | GGAGCCACTATCGACTAC | 50 |
| TAGC1029 | pBR3 | TCCCCATCGGTGATGTC | 50 |
| TAGC1030 | PBR4 | CTCGGCACCGTCACC | 50 |
| TAGC1031 | pc3.1GFP-Topo-RP | CCCATTAACATCACC | 50 |
| TAGC1032 | pcDNA1.1-RP | CTCTGTAGGTAGTTTGTC | 50 |
| TAGC1033 | pcDNA1.1-RP/1 | TAGAATCAGTAGTTTAACAC | 50 |
| TAGC1034 | pcDNA3.1-FP | CTCTGGCTAACTAGAGAAC | 50 |

| | | | |
|----------|---------------|----------------------------|----|
| TAGC1035 | pcDNA3.1-RP/1 | CAAACAACAGATGGCTGGC | 50 |
| TAGC1036 | PCR1 | CGGGCCTCTTCGCTATT | 50 |
| TAGC1037 | PCR1minus1 | GGGCCTCTTCGCTATT | 50 |
| TAGC1038 | PCR2 | TTAGCTCACTCATTAGG | 50 |
| TAGC1039 | pDONOR-FP | TAACGCTAGCATGGATCTC | 50 |
| TAGC1040 | pDONOR-RP | GCAATGTAACATCAGAGAT | 50 |
| TAGC1041 | pEGFP-FP | TTTAGTGAACCGTCAGATC | 50 |
| TAGC1042 | pEGFP-N | CCGTCCAGCTCGACCAG | 50 |
| TAGC1043 | pEGFP-RP | AACAGCTCCTCGCCCTTG | 50 |
| TAGC1044 | pEGFP_C2-FP | GATCACATGGTCCTGCTG | 50 |
| TAGC1045 | pEGFP_C2-RP | TTTAAAGCAAGTAAAACCTC | 50 |
| TAGC1046 | pESP-RP | TCCAAAAGAAGTCGAGTGG | 50 |
| TAGC1047 | pET-24a | GGGTTATGCTAGTTATTGCTCAG | 50 |
| TAGC1048 | pET-RP | CTAGTTATTGCTCAGCGG | 50 |
| TAGC1049 | pET30er-FP | TCATCATTCTTCTGGTCTG | 50 |
| TAGC1050 | pET36b-FP | GCTAAATTCGAACGCCAG | 50 |
| TAGC1051 | pETBlue-RP | AATAGCTTTAATGCGGTAG | 50 |
| TAGC1052 | pGBT9-FP | AGTGCGACATCATCATCG | 50 |
| TAGC1053 | pGBT9-RP | CGTTTTAAAACCTAAGAGTCAC | 50 |
| TAGC1054 | pGEX-3 | GGAGCTGCATGTGTCAGAG | 50 |
| TAGC1055 | pGEX-5 | CTGGCAAGCCACGTTTGG | 50 |
| TAGC1056 | pGEX3-RP | TCAAGAATTATACACTCCG | 50 |
| TAGC1057 | pGEX5-FP | AACGTATTGAAGCTATCCC | 50 |
| TAGC1058 | pGFP-FP | GAGTTTGTAACAGCTGCTG | 50 |
| TAGC1059 | pGL2 | TCTTTATGTTTTTGGCGTC | 50 |
| TAGC1060 | pGL3-1 | GAGCTGACTGGGTTGAAG | 50 |
| TAGC1061 | pGL3-RV3 | CTAGCAAATAGGCTGTCC | 50 |
| TAGC1062 | pIRES-RP | TATAGACAAACGCACACCG | 50 |
| TAGC1063 | pJet1-FP | ACTACTCGATGAGTTTTTCGG | 50 |
| TAGC1064 | pJet1-RP | TGAGGTGGTTAGCATAGTTC | 50 |
| TAGC1065 | pMalE | TCAGACTGTTCGATGAAGC | 50 |
| TAGC1066 | pMT-FP | CATCTCAGTGCAACTAAAG | 50 |
| TAGC1067 | PolyA-B | AAAAAAAAAAAAAAAAAAAAAAAAAB | 50 |
| TAGC1068 | PolyC-D | CCCCCCCCCCCC | 50 |
| TAGC1069 | PolyT-V | TTTTTTTTTTTTTTTTTTTTTV | 50 |
| TAGC1070 | pQE-FP | CGGATAACAATTTACACAG | 50 |
| TAGC1071 | pQE-RP | GTTCTGAGGTCATTACTGG | 50 |
| TAGC1072 | pREP-FP | GCTCGATACAATAAACGCC | 50 |
| TAGC1073 | pRH-FP | CTGTCTCTATACTCCCCTATAG | 50 |
| TAGC1074 | pRH-RP | CAAAATTCAATAGTTACTATCGC | 50 |
| TAGC1075 | pRSET-RP | ATGCTAGTTATTGCTCAGC | 50 |
| TAGC1076 | pRSET-Rpnew | GGGTTATGCTAGTTATTGC | 50 |
| TAGC1077 | pTal-FP | CGGGAGGTACTTGGAGCG | 50 |
| TAGC1078 | pTeSp-1 | CCTCCATAGAAGACACC | 50 |
| TAGC1079 | pTeSp-2 | GCCCTGCCACTCATCG | 50 |
| TAGC1080 | pTI2-1-1 | GGCAAATATTCTGAAATGAGC | 50 |
| TAGC1081 | pTI2-1-2 | GCGTTTCACTTCTGAGTTCG | 50 |
| TAGC1082 | pTrcHis-FP | GACCGGAATTATCGATTAAC | 50 |
| TAGC1083 | pTrcHis-RP | CTGATTTAATCTGTATCAGG | 50 |
| TAGC1084 | pTYB-FP | GGAAGACGATTATTATGGG | 50 |
| TAGC1085 | pUC-F | GCCAGTGAATTCGAGCTCGG | 50 |

| | | | |
|----------|-------------|---------------------------|----|
| TAGC1086 | pUC-R | TGCCTGCAGGTCGACTCTAG | 50 |
| TAGC1087 | pVL1 | GGGTTTAACATTACGGATTTTC | 50 |
| TAGC1088 | pVL2 | GACGCACAACTAATATCAC | 50 |
| TAGC1089 | pVP22-FP | CGCTTCTCGCCCCAGAC | 50 |
| TAGC1090 | RPC | ACAGCTATGACCATGATTACG | 50 |
| TAGC1091 | SK | CGCTCTAGAAGTACTAGTGGATC | 50 |
| TAGC1092 | SP6 | ATTTAGGTGACACTATAGAA | 50 |
| TAGC1093 | T3 | ATTAACCCTCACTAAAGGGA | 50 |
| TAGC1094 | T7 | TAATACGACTCACTATAGGG | 50 |
| TAGC1095 | T7-981079 | TAATACGACTCACTATAG | 50 |
| TAGC1096 | T7minus1 | AATACGACTCACTATAGGG | 50 |
| TAGC1097 | TET-1 FP-1 | GGGTGCGCATGATCCTCTAGAGT | 50 |
| TAGC1098 | TET-1 RP-1 | TAAATTGCACTGAAATCTAGAAATA | 50 |
| TAGC1099 | Topo-1 | TCGGATCCACTAGTAACG | 50 |
| TAGC1100 | Topo-2 | GTGTGATGGATATCTGC | 50 |
| TAGC1101 | Topo-2N | AGATGCATGCTCGAGCGG | 50 |
| TAGC1102 | TriplEx5`LD | GAAGCGCGCCATTGTGTTGG | 50 |
| TAGC1103 | UP-2 | CGTTGTAAAACGACGGCC | 50 |
| TAGC1104 | UP-40 | GTTTTCCAGTCACGAC | 50 |

Price list valid from 19th Sept. 2012

Prices are in DKK excl. VAT

Catalog prices are subject to change without notice.

How to order:

Please contact our customer service when ordering standard primers.